



Luxator® all you need for a successful extraction!

The ability to reduce trauma during tooth extraction has become extremely important, especially for patients requiring subsequent implant placement with minimal bone lost. Modern implant techniques require real alternatives to traumatic forcep extractions.

The Luxator product line offers a safe and precise solution for successful extractions. The fine tapered blade compresses the alveolar, cuts the membrane, and gently eases the tooth from its socket reducing damage to surrounding tissue, keeping a better anatomy for an implant site. The shape and size of the handle is designed to minimize the force needed during the extraction. The thin blade can be easily inserted between the bone and the root with minimal destruction.

Manufactured with Swedish stainless steel, these instruments are extremely durable, retaining their strength and sharpness after repeated use and sterilization. With their acclaimed ergonomic handle Luxator Instruments are discernibly different and recognized worldwide as one of the best surgical instruments on the market!

Luxator® is a registered trademark owned by Directa.



Luxator®

Luxator® Periotomes page: 4-29

page: 8-11 page: 12-15 page: 16-19 page: 20-23 page: 24-29 Luxator Periotome Luxator Short Periotome Luxator Dual Edge Periotome Luxator Titanium Periotome Luxator LX Mechanical Periotome

Luxator® Forte Elevator page: 30-33

Luxator® RootPicker page: 36-39



Preserved bone integrity and dramatically less trauma for your patients and their dentition Luxator Periotomes are specially designed with fine tapering blades that compress the alveolar bone, cut the membrane and gently ease the tooth from its socket. The whole operation is performed with a minimum of tissue damage.





Ergonomic design, world acclaimed handle design

The Luxator instruments were invented by a Swedish dentist to make extractions as trauma free as possible. He developed subtleties in the design that only a practicing dentist would appreciate. With an acclaimed, ergonomic handle Luxator Instruments are discernibly different. The Luxator is ergonomically designed for comfort and control during use. The shape and size of the handle is designed to minimize the force needed. The Luxator extraction instrument product range enhances the outstanding selection of high quality extraction instruments from Directa.

Reasons to use Luxator® Periotome

Ergonomic design

World acclaimed ergonomic design of the handle. Minimizes fatigue.

Atraumatic extraction

The fine tapered blade compresses the alveolar bone and cuts the periodontal ligaments.

Maximum tactility

The use of highest quality materials gives the instruments outstanding tactility.

Wide range of sizes and shapes

Unique variety of sizes and shapes enables users to find the ideal instrument for any extraction.





periodontal ligament knife

Luxator Periotomes are the preferred instruments for performing extractions. Luxator Periotomes are specially designed periodontal ligament knives with fine tapering blades that compress the alveolar bone, cut the ligament and gently ease the tooth from its socket. The whole operation

is performed with a minimum of tissue damage. Available in 8 blade sizes. Autoclavable Fixus 1

Prod No: 506352, L2S, 2mm, straight blade, dark green Prod No: 506330, Luxator Kit K4 Prod No: 506331, Luxator Kit K7 Prod No: 506332, Luxator Kit K VET

case included.







L2S
Periotome
506352
2 mm
Straight
Standard



L3S Periotome 506340 3 mm Straight Standard



L3C Periotome 506341 3 mm Curved Standard







L3IC Periotome 506354 3 mm Inverted Curved Standard



L5S Periotome 506342 5 mm Straight Standard



L5C Periotome 506342 5 mm Curved Standard

Reasons to use Luxator® Short Periotome

Shorter blade

World acclaimed handle design. Minimizes fatigue. The 10mm shorter blade enables an optimal working position for smaller hands.

Atraumatic extractions

The fine tapered blade compresses the alveolar bone and cuts the periodontal ligaments.

Maximum tactility

The use of highest quality material gives the instruments outstanding tactility.

Wide range of sizes and shapes

Five sizes and shapes enable users to find the ideal instrument for any extraction.





Luxator® Short Periotome the original with a shorter blade

Prod No: 506361, S3CA, 3mm, contra angle, dark blue Prod No: 506333, Luxator Kit SK4

Attains the correct working position - even with smaller hands. The shaft is 10mm shorter than on regular Luxator Periotome. Preserved bone integrity and dramatically less trauma for your patients and their dentition. Available in 5 blade sizes. Autoclavable Fixus 1 case included.









S2S Periotome 506358 2 mm Straight Short

S3S Periotome 506360 3mm Straight Short

S3C Periotome 506359 3 mm Curved Short

S3CA
Periotome
506361
3 mm
Contra Angle
Short



\$5\$
Periotome
506362
5 mm
Straight
Short

Reasons to use **Luxator® Dual Edge**

Dual EdgeIntegrated periotome and Luxator tip.

Safer insertion

The Dual Edge tip offers safer stepped insertion and reduces the risk of slipping.

Ideal for fractured teeth

luxating procedure.

Perfect starter in difficult ases

The thin periotome part is ideal to use as a starter for difficult extractions.





Luxator® Dual Edge Periotome innovative dual purpose tip

Prod No: 506356 Luxator DE3, 3mm/1,5mm, dual edge, tan

Integrated probing Luxator Periotome. Safer stepped insertion reduces the risk of slipping. Ideal for deeply fractured teeth and decayed roots. Available in 2 blade sizes. Autoclavable Fixus 1 case included.



DE3 Dual Edge 506356 3/1.5 mm Straight Standard



DE5 Dual Edge 506357 5/3 mm Straight Standard



Reasons to use Luxator® Titanium

Stays sharpThe durable titanium coating ensures that the intsument stays sharp.

Low friction

The titanuim coating lowers the friction with 43% compared with stainless steel, thus esier to insert into the socket.

Ergonomic design

World acclaimed ergonomic design of the handle minimizes fatigue.

Wide range of sizes and shapes Unique variety of sizes and shapes enables users to find the ideal instrument for any extraction.





Luxator® Titanium Periotome sharpen-free instrument

Prod No: 506433 Luxator L3S TiN, grey

Luxator Periotomes are specially designed with fine tapering blades that compress the alveolar bone, cut the ligament and gently ease the tooth from its socket.

The Luxator Titanium Periotome offers enhanced durability and extractions are performed with a minimum of tissue damage. The titanium coating lowers the friction with 43% compared with stainless steel, making it easier to glide into the socket. The new durable titanium coated tip does not require sharpening. Available in 4 blade sizes. Autoclavable Fixus 1 case included.







L3S TiN Periotome 506433 3mm Straight Titanium



DE3 TiN Dual Edge 506435 3 mm Straight Titanium



L3A TiN Periotome 506436 3 mm Angled Titanium



Reasons to use Luxator® LX

Optimal access

Undisputable access in hard to reach areas.

Atraumatic extractions

Minimal risk of tissue damage. Fast recovery for the patient. Facilitates direct placement of

Vertical reciprocating tip Enables safe cutting of the periodontal ligaments.

Self directing tip

Allows the periotome to follow the root surface and reduce the risk of












Luxator® LX mechanical periotome

Prod No: 506443 Luxator LX Starter Kit

Setting a new standard for trauma free extraction, the Luxator LX mechanical periotome offers improved safety, precision and access.

Safer extractions with minimal tissue damage, the Luxator LX Mechanical Peritome helps access hard to reach areas. The self directing, reciprocating movement allows the periotome tip to follow the root surface and enables a secure cutting of the fibers in a patient friendly way, reducing the risk of bone damage. With a reciprocating tip The Luxator LX Mechanical Periotome will only cut when under pressure, which means no risk of slipping. The titanium coated tip is a durable tip and remains sharp.



Luxator® LX periotome tips

LX2S Short length 15 mm, width 2 mm length 15 mm width 3 mm

LX2L Long length 19 mm width 2 mm length 19 mm width 3 mm

The Periotome tips are made of titanium coated stainless steel. Do not sharpen.



Luxator® LX hand piece

- Standard E-fitting
 Anodized aluminium
 Recommended speed 1000-4000 rpm (low speed)
 Maximum speed 4000 rpm
 Axial stroke

Caution

Do not exceed 4000 rpm. Exceeding could cause discomfort in terms of vibrations to patient or operator.



Insert the Periotome Tip into the handpiece, starting with 2mm tip.



Activate the mechanical periotome and luxate at a low speed (2000 - 4000 rpm). $\,$



Let the Periotome Tip work around the tooth to cut the periodontal fibers. Press firmly and follow the surface of the root.



Use a Luxator Forte to judge the actual luxation status. If necessary continue with Luxator LX. If entering is hard, try to get an access point with a thin Luxator Periotome and then continue the luxation with Luxator LX.



Remove the tooth with ordinary forceps or Luxator RootPicker.



Successful extraction with minimal tissue damage.



Luxator® Forte Elevator for controlled and safer tooth elevation

Luxator Forte Elevator incorporates the world acclaimed ergonomic handle and a unique cornered, rigid tip that enables elevation of roots without risk of damaging the instrument. Luxator Forte is the perfect match to Luxator Periotome to end a successful extraction!





Reasons to use Luxator® Forte

Ergonomic design World acclaimed handle. Ergonomic

design. Minimizes fatigue.

Thin,
yet strong
The use of highest quality materials allows the tip to be very thin without its strenght thus offering better access.

Cornered

tip designThe unique cornered design grasps the tooth for safe and controlled elevations.

Maximum tactility

Highest quality materials for instruments with outstanding tactility.





Incorporates the world acclaimed ergonomic handle and a cornered, rigid tip that enables elevation of roots without risk of damaging the instrument. The unique cornered tip design grasps dentine or root cementum for controlled and safer tooth elevation. Available in 4 blade sizes. Autoclavable Fixus 1 case included.









F25 Forte 506370 2.5 mm Straight Standard

F32 Forte 506371 3.25 mm Straight Standard

F32C Forte 506373 3,25 mm Curved Standard

F40 Forte 506372 4 mm Straight Standard



Reasons to use Luxator® RootPicker

Rigid and slim

Designed to offer a strong grip with a slim design that is easy to insert into the alveolus.

Unique design
Ideal for handling components used in surgery, such as implant accessories.

Four Corner gripProvides an excellent grip on root tips.

Straight or angled

For best access.





grips debris after extraction

Prod No: 506450 RootPicker Straight Prod No: 506451 RootPicker Angled Prod No: 506452 RootPicker Assortment

During extractions annoying tooth fragments or loose root tips remain attached to the periodontal ligament. These pieces are very small and slippery, making it difficult to grip them with ordinary tweezers or forceps. RootPicker instruments address these problems with the necessary properties to complete the extraction successfully. Luxator RootPicker is a short and rigid instrument easy to insert into the alveolus. The cleverly designed tip provides a sure grip on the debris. Their straight and angled tips make them ideal for use in the anterior area as well as in the posterior maxillary and mandibular areas.



Luxator® RootPicker grips implantology components

RootPicker also provides excellent grip for the various components used in surgery such as healing abutments, parallelism indicators, screwdrivers and posts for immediate prosthesis. Their straight and angled design make them ideal for use in the anterior area as well as in the posterior maxillary and mandibular areas.





Design By Dentists™

Directa's products are developed and evaluated by our selected team of highlyqualified and renowned Swedish and international dental professionals with the aim of making life easier for the dental professionals in the daily clinical work.

DIRECTA directadental.com

DIRECTA AB

P.O. Box 723, 194 27 Upplands Väsby, SWEDEN +46 8 506 505 75

DIRECTA INC

64 Barnabas Road, Unit 3 Newtown, CT 06470, United States 203.491.2273

DIRECTA GERMANY

Hammergutstr. 11 92245 Kümmersbruck, GERMANY

info@directadental.com www.directadental.com www.directadentalgroup.com www.facebook.com/directadentalgroup www.instagram.com/directadentalgroup www.twitter.com/D_DentalGroup www.youtube.com/user/DirectaAB

DirectaDentalGroup

DIRECTA topdental orsing parkell

