



W o r l d L e a d e r i n S p h e r i c a l A t t a c h m e n t s

# CODE LIST

for Dentist and Dental Laboratory



2017

w w w . r h e i n 8 3 3 . c o m



## CASTABLE



Ref. 092ECQ

### OT EQUATOR CASTABLE

**KIT CONTAINS:**

- 2 SINGLE MALES
- 2 STAINLESS STEEL HOUSINGS
- 4 RETENTIVE CAPS  
(2 SOFT, 2 STANDARD)

## IMPLANTS



Ref. 130

### OT EQUATOR FOR IMPLANTS

**KIT CONTAINS:**

- 1 TITANIUM ABUTMENT
- 1 TITANIUM HOUSING
- 1 PROTECTIVE DISK
- 4 RETENTIVE CAPS  
(1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)



Ref. 030

- 1 TITANIUM ABUTMENT

## BARS



Ref. 160EQB

### OT EQUATOR THREADED SLEEVE SYSTEM FOR BONDING

**KIT CONTAINS:**

- 2 EQUATOR THREADED for titanium sleeve (1.6 mm THREAD)
- 2 THREADED SLEEVES (1.6 mm THREAD)
- 2 STAINLESS STEEL HOUSINGS
- 2 EQUATOR SPACERS
- 8 RETENTIVE CAPS  
(2 EXTRA-SOFT, 2 SOFT, 2 STANDARD, 2 PROCESSING)

## ELASTIC SEEGER PASSIVE BAR CONNECTIONS SYSTEM



Ref. 158ESA

### ELASTIC SEEGER

**KIT CONTAINS:**

- 1 CASTABLE CYLINDER FOR SEEGER
- 1 SELF-EXTRACTING SEEGER
- 1 TITANIUM SCREW FOR SELF-EXTRACTING SEEGER

## CAPS ASSORTMENT KIT



Ref. 192ECE

### CAPS ASSORTMENT KIT

**KIT CONTAINS:**

- 1 TITANIUM HOUSING
- 1 PROCESSING CAP
- 1 PROTECTIVE DISK
- 4 RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)

## SMART BOX KIT



Ref. 335SBC

### SMART BOX KIT

**KIT CONTAINS:**

- 1 SMART BOX HOUSING WITH BLACK CAP FOR LABORATORY
- 1 PINK PROTECTIVE DISK
- 4 RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)

## SPARE PARTS



|                     |   |
|---------------------|---|
| <b>Ref. 151SS</b>   | • 2 CASTABLE MALES                                    |
| <b>Ref. 140CEV</b>  | • 4 RETENTIVE CAPS - VIOLET "STRONG" 2700g            |
| <b>Ref. 140CET</b>  | • 4 RETENTIVE CAPS - CLEAR "STANDARD" 1800g           |
| <b>Ref. 140CER</b>  | • 4 RETENTIVE CAPS - PINK "SOFT" 1200g                |
| <b>Ref. 140CEG</b>  | • 4 RETENTIVE CAPS - YELLOW "EXTRA-SOFT" 600g         |
| <b>Ref. 140CEN</b>  | • 4 PROCESSING CAPS - BLACK (FOR LABORATORY USE ONLY) |
| <b>Ref. 141CAE</b>  | • 2 STAINLESS STEEL HOUSINGS                          |
| <b>Ref. 141CTE</b>  | • 2 TITANIUM HOUSINGS                                 |
| <b>Ref. 330SBE</b>  | • 1 SMARTBOX HOUSING WITH BLACK CAP                   |
| <b>Ref. 335CSB</b>  | • 4 BLACK SMARTBOX CAPS (FOR LABORATORY USE ONLY)     |
| <b>Ref. 039SFE2</b> | • 1 EQUATOR THREADED 2.0mm (DIRECT INTO 2mm THREAD)   |
| <b>Ref. 339SFE</b>  | • 1 EQUATOR THREADED 1.6mm (FOR TITANIUM SLEEVE)      |
| <b>Ref. 239GSF</b>  | • 1 TITANIUM SLEEVE (1.6mm THREAD)                    |
| <b>Ref. 339SSE</b>  | • 1 OT EQUATOR SPACER                                 |
| <b>Ref. 158VAE</b>  | • 1 TITANIUM SCREW FOR SELF-EXTRACTING SEEGER         |
| <b>Ref. 158SAE</b>  | • 1 SELF-EXTRACTING WHITE SEEGER (FOR BAR LOCKING)    |
| <b>Ref. 158CCE</b>  | • 1 CASTABLE CYLINDER height 2,5mm (FOR SEEGER)       |
| <b>Ref. 15CCA</b>   | • 1 HIGH CASTABLE CYLINDER height 3,5mm (FOR SEEGER)  |

## ACCESSORIES



|                     |   |
|---------------------|---|
| <b>Ref. 144TTE</b>  | • 1 OT EQUATOR TITANIUM TRANSFER + TITANIUM SCREW   |
| <b>Ref. 144MTE</b>  | • 2 IMPRESSION COPINGS (FOR INDIVIDUAL TRAY)  |
| <b>Ref. 044CAIN</b> | • 2 IMPRESSION COPINGS NORMAL SIZE - OT EQUATOR (WITH INTERCHANGEABLE CAP FOR PICK UP IMPRESSION) |
| <b>Ref. 144AE</b>   | • 2 LABORATORY ANALOGS  |

## TOOLS



|                    |   |
|--------------------|---|
| <b>Ref. 74AC01</b> | • 1 PARALLELOMETER KEY "NORMAL SIZE"  |
| <b>Ref. 774CHE</b> | • 1 OT EQUATOR SQUARE SCREW DRIVER + HOLDER (SQUARE 1.25mm)   |
| <b>Ref. 774HC</b>  | • 1 OT EQUATOR INTERCHANGEABLE HOLDER   |
| <b>Ref. 760CE</b>  | • 1 SQUARE DRIVER CONNECTOR for torque CONTROLLER CONTRA ANGLE (SQUARE 1.25mm)  |
| <b>Ref. 485IAC</b> | • 1 MULTIUSE CAPS INSERTER/EXTRACTOR TOOL   |
| <b>Ref. 185SIS</b> | • 1 METAL INSERTION TOOL FOR SEEGER USE TOGETHER WITH 124ICP BLU HANDLE   |
| <b>Ref. 491EC</b>  | • 1 MULTIUSE EXTRACTOR TOOL FOR CAPS (WITH HOLE STEM HOLDER)  |
| <b>Ref. 760CRD</b> | • 1 RATCHET TORQUE CONTROL DEVICE (FOR OT EQUATOR, SPHERO BLOCK AND FLEX 15/35NCM STRENGHT - MAX 50NCM TORQUE, SUGGESTED 25NCM) |

## OT UNILATERAL

EXTRA CORONAL ATTACHMENTS WITH COMBINED FUNCTION



Ref. 098UNL

### OT UNILATERAL + CONNECTOR

#### KIT CONTAINS:

- 1 UNILATERAL CASTABLE ATTACHMENT
- 1 UNI-BOX CASTABLE
- 1 PLASTIC POSITIONER MICRO SIZE
- 1 PINK CAP MICRO SIZE
- 1 STRATEGY YELLOW CAP (for duplication technique)
- 1 BLACK CAP MICRO SIZE (for laboratory use only)
- 1 STRATEGY BLACK CAP (for laboratory use only)
- 1 CONNECTOR



Ref. 156BCU

- 1 UNILATERAL CASTABLE ATTACHMENT



Ref. 156CCU

- 1 UNI-BOX CASTABLE

## OT CAP TECNO



NORMAL

Ref. 093CTN

### OT CAP TECNO NORMAL/MICRO

#### KIT CONTAINS:

- 2 PRE-ANGULATED CASTABLE EXTENSIONS
- 2 TITANIUM SINGLE THREADED SPHERES
- 2 TITANIUM BOND-IN THREADED SLEEVES
- 2 TRANSPARENT CAPS NORMAL/MICRO
- 2 PINK CAPS NORMAL/MICRO
- 2 YELLOW CAPS NORMAL/MICRO



MICRO

Ref. 093CTM



Ref. 193BCA

- 2 PRE-ANGULATED CASTABLE EXTENSIONS

## SINGLE THREADED SPHERES

WITH THREADED SLEEVE TO BE BONDED NORMAL - MICRO



Ref. 139KSFN

### KIT SINGLE THREADED SPHERE NORMAL SIZE + TiN WITH SLEEVE TO BE BONDED

#### KIT CONTAINS

- 2 TITANIUM NORMAL SINGLE THREADED SPHERE  
(*HEX mm. 1.3 THREAD mm. 1.6*)
- 2 TITANIUM BOND-IN THREADED SLEEVE
- 2 S/STEEL SPACER for the THREADED SLEEVE (Normal sphere)



Ref. 139KSFM

### KIT SINGLE THREADED SPHERE MICRO SIZE + TiN WITH SLEEVE TO BE BONDED

#### KIT CONTAINS

- 2 TITANIUM MICRO SINGLE THREADED SPHERE  
(*HEX mm. 0.9 THREAD mm. 1.6*)
- 2 TITANIUM BOND-IN THREADED SLEEVE
- 2 S/STEEL SPACER for the THREADED SLEEVE (Micro sphere)

## SPARE PARTS



Ref. 239SFN

- 1 SINGLE THREADED SPHERE NORMAL SIZE + TiN  
FOR THREADED SLEEVE HEX mm. 1.3 THREAD mm. 1.6



Ref. 239SFM

- 1 SINGLE THREADED SPHERE MICRO SIZE + TiN  
FOR THREADED SLEEVE HEX mm. 0.9 THREAD mm. 1.6



Ref. 239GSF

- 1 TITANIUM BOND-IN THREADED SLEEVE  
FOR NORMAL AND MICRO - mm. 1.6



Ref. 239SSN

- 1 **STAINLESS STEEL SPACER for THREADED SLEEVE (Normal sphere)**

Ref. 239SSM

- 1 **STAINLESS STEEL SPACER for THREADED SLEEVE (Micro sphere)**

## SINGLE THREADED SPHERES

FOR CAD/CAM OR SELF TAPPING 2.0mm THREAD - NORMAL AND MICRO SIZE



Ref. 039SFN2

- 1 **SINGLE THREADED SPHERE NORMAL SIZE + TiN**  
HEX 1.3mm, THREAD 2.0mm

Ref. 039FM2

- 1 **SINGLE THREADED SPHERE MICRO SIZE + TiN**  
HEX 0.9mm, THREAD 2.0mm



Ref. 760CBR

- 1 **HEX DRIVER FOR THREADED SPHERE NORMAL**  
HEX 1.3mm



Ref. 772CSF

- 1 **HEX DRIVER FOR THREADED SPHERE MICRO**  
HEX 0.9mm

## OT CAP IMPRESSION COPINGS



Ref. 044CAIN

- 2 **IMPRESSION COPINGS NORMAL SIZE**  
Ø2,5mm WITH INTERCHANGEABLE CAP



Ref. 044CAI22

- 2 **IMPRESSION COPINGS**  
Ø1,25mm SPHERES WITH INTERCHANGEABLE CAP



Ref. 044CAIM

- 2 **IMPRESSION COPINGS MICRO SIZE**  
Ø1,8mm WITH INTERCHANGEABLE CAP

## ATTACHMENTS OT CAP & OT BOX - OT STRATEGY - OT BAR OT VERTICAL - OT TECNO - TITANIUM SINGLE SPHERE

MASTER STARTER KIT FOR LABORATORY



Ref. 005SKLUS

### KIT CONTAINS:

#### OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (normal and micro size)
- 4 castable spheres (normal and micro size)
- 2 titanium single spheres (normal and micro size) + sphere holder + 2 sliding base
- 4 castable OT CAP bars (normal and micro size)
- 2 complete castable OT CAP TECNO bars (normal and micro size)
- 4 castable OT BOX CLASSIC bars (upper + lower part) normal and micro size
- 2 castable OT BOX SPECIAL bars (normal and micro size)
- 6 castable OT BOX connectors
- 4 castable OT BOX MONO housings (normal and micro size)
- 4 plastic positioner rings (normal and micro size)
- 16 OT CAP retentive caps (normal and micro size)  
Clear - standard retention - PINK soft retention - YELLOW extra soft retention  
GREEN very elastic retention

- 4 TITAN CAP titanium caps (normal and micro size)
- 4 BLACK CAPS- only for lab use - (normal and micro size)
- 8 S/STEEL housings
- 4 ANALOGUES for plaster models (normal and micro size)
- 4 protective disks
- 4 aluminium spacers (normal and micro size)

#### OT STRATEGY:

- 4 castable spheres (2 standard and 2 high)
- 2 castable STEADY
- 6 retentive caps: 2 clear - 2 pink - 2 yellow for housing
- 2 black caps for housing only for lab use
- 2 s/steel housings
- 2 plastic positioner rings
- 6 retentive caps for duplication technique:  
Clear - standard retention - PINK soft retention - YELLOW extra soft retention
- 2 black caps for duplication technique only for lab use
- 2 ANALOGUES for plaster models

#### OT BAR MULTIUSE:

- 2 castable bar
- 1 gingival connector
- 2 positioner clip A
- 2 positioner clip B
- 4 castable housings
- 4 clips - 2 PINK - soft retention - 2 YELLOW medium retention

#### OT VERTICAL:

- 2 castable male
- 2 castable STEADY
- 4 clips - 2 WHITE- standard retention - 2 GREEN soft retention
- 2 castable parallelometer keys
- 2 ceramic bars

#### TOOLS:

- 1 Pincers
- 2 Reamer tools for caps (normal and micro size)
- 1 Reamer tools for OT STRATEGY caps
- 2 Parallelometer keys for OT CAP (normal and micro size)
- 1 Parallelometer key for OT STRATEGY size
- 2 Parallelometer keys for OT CAP TECNO (normal and micro size)
- 1 Parallelometer key for OT BAR MULTIUSE
- 1 Blue Plastic "multiuse" insertion tool  
(for steel tool and for inserting the caps in the mouth)
- 1 OT CAP insertion tool
- 1 OT STRATEGY insertion tool
- 1 OT BAR clips insertion tool
- 1 OT VERTICAL clips insertion tool

# ATTACHMENTS OT CAP & OT BOX - OT STRATEGY - OT BAR

## OT VERTICAL - OT UNILATERAL - OT EQUATOR

BASIC LABORATORY KIT



Ref. 005SKLBUS

### KIT CONTAINS:

#### OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (normal and micro size)
- 8 castable spheres normal/micro
- 2 castable OT CAP bars (normal and micro size)
- 2 castable OT BOX CLASSIC bars (upper + lower part) normal and micro size
- 2 castable OT BOX SPECIAL bars (normal and micro size)
- 6 castable OT BOX connectors
- 4 castable OT BOX MONO housings (normal and micro size)
- 8 plastic positioner rings (normal and micro size)
- 13 OT CAP assorted retentive caps (normal and micro size)  
PINK soft retention - YELLOW extra soft retention
- 13 BLACK CAPS- only for lab use - (normal and micro size)
- 4 S/STEEL housings (normal and micro size)

#### OT STRATEGY:

- 4 castable spheres (2 standard and 2 high)
- 2 castable STEADY
- 2 pink retentive caps
- 2 s/steel housings
- 2 plastic positioner rings
- 3 retentive caps for duplication technique:  
PINK soft retention - YELLOW extra soft retention
- 4 black caps for duplication technique only for lab use

#### OT EQUATOR:

- 2 single males
- 2 stainless steel housings
- 4 retentive caps (2 soft, 2 standard)
- 2 processing black caps

#### OT UNILATERAL:

- 1 unilateral castable attachment
- 1 uni-box castable
- 1 plastic positioner micro size
- 2 ot cap micro size: 1 pink - 1 black processing cap
- 2 ot strategy (for duplication technique): 1 yellow - 1 black processing cap
- 1 connector

#### OT BAR MULTIUSE:

- 1 castable bar
- 1 gingival connector
- 4 positioner clip (A - B)
- 2 castable housings
- 4 clips - 2 PINK - soft retention - 2 YELLOW medium retention

#### OT VERTICAL:

- 2 castable male
- 2 castable STEADY
- 4 clips - 2 WHITE- standard retention - 2 GREEN soft retention
- 2 castable parallelometer keys
- 2 ceramic bars

#### TOOLS:

- 1 Tweezers
- 1 Parallelometer keys for OT CAP normal
- 1 Parallelometer keys for OT CAP micro
- 1 Parallelometer key for OT STRATEGY
- 1 Parallelometer key for OT BAR MULTIUSE
- 1 Blue Plastic "multiuse" insertion tool  
(for steel tool and for inserting the caps in the mouth)
- 1 OT CAP insertion tool
- 1 OT STRATEGY insertion tool
- 1 OT BAR clips insertion tool
- 1 OT VERTICAL clips insertion tool

## OT CAP SYSTEM ATTACHMENTS




Ref. 038STF  
NORMAL SIZE

### TITANIUM FLEX SINGLE SPHERE NORMAL SIZE + TiN FOR LASER SOLDERING - FOR BONDING

#### KIT CONTAINS:

- 2 SINGLE TITANIUM FLEX SPHERES
- 2 PINK CAPS - SOFT RETENTION
- 1 SPHERE HOLDER
- 2 CASTABLE SLIDING BASES



Ref. 038STN  
NORMAL SIZE

### TITANIUM SINGLE SPHERE NORMAL SIZE + TiN FOR LASER SOLDERING - FOR BONDING

#### KIT CONTAINS:

- 2 SINGLE TITANIUM SPHERES
- 2 PINK CAPS - SOFT RETENTION
- 1 SPHERE HOLDER
- 2 CASTABLE SLIDING BASES



Ref. 038STM  
MICRO SIZE

### TITANIUM SINGLE SPHERE MICRO SIZE + TiN FOR LASER SOLDERING - FOR BONDING

#### KIT CONTAINS:

- 2 SINGLE TITANIUM SPHERES
- 2 PINK CAPS - SOFT RETENTION
- 1 SPHERE HOLDER
- 2 CASTABLE SLIDING BASES



Ref. 092CAN

### NORMAL SIZE OT CAP KIT CONTAINS

- 4 CASTABLE SINGLE SPHERES
- 4 PINK RETENTIVE CAPS
- 4 STAINLESS STEEL HOUSINGS
- 4 PLASTIC POSITIONERS



Ref. 092CAM

### MICRO SIZE OT CAP KIT CONTAINS

- 4 CASTABLE SINGLE SPHERES
- 4 PINK RETENTIVE CAPS
- 4 STAINLESS STEEL HOUSINGS
- 4 PLASTIC POSITIONERS



Ref. 196BCN

### "ECONOMIC" OT CAP NORMAL SIZE KIT CONTAINS

- 1 CASTABLE PRESECTIONED BAR
- 1 CASTABLE PRESECTIONED BEVELED BAR
- 4 CLEAR RETENTIVE CAPS (Standard retention)



Ref 197BCM

### "ECONOMIC" OT CAP MICRO SIZE KIT CONTAINS

- 1 CASTABLE PRESECTIONED BAR
- 1 CASTABLE PRESECTIONED BEVELED BAR
- 4 CLEAR RETENTIVE CAPS (Standard retention)



Ref. 099BSN

### OT CAP & Mono OT BOX FOR FRAME NORMAL SIZE KIT CONTAINS

- 2 CASTABLE PRESECTIONED BARS (1 straight + 1 beveled)
- 4 CLEAR RETENTIVE CAPS
- 4 CASTABLE MONO OT BOX
- 4 PLASTIC POSITIONERS



Ref. 099BSM

### OT CAP & Mono OT BOX FOR FRAME MICRO SIZE KIT CONTAINS

- 2 CASTABLE PRESECTIONED BARS (1 straight + 1 beveled)
- 4 CLEAR RETENTIVE CAPS
- 4 CASTABLE MONO OT BOX
- 4 PLASTIC POSITIONERS



Ref. 058BSN

### OT BOX SPECIAL NORMAL SIZE + CONNECTORS KIT CONTAINS

- 2 OT BOX SPECIAL BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS



Ref. 058BSM

### OT BOX SPECIAL MICRO SIZE + CONNECTORS KIT CONTAINS

- 2 OT BOX SPECIAL BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS



Ref. 153BCN

### OT BOX CLASSIC NORMAL SIZE + CONNECTORS KIT CONTAINS:

- 2 UPPER BARS
- 2 LOWER BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS



Ref. 153BCM

### OT BOX CLASSIC MICRO SIZE + CONNECTORS KIT CONTAINS:

- 2 UPPER BARS
- 2 LOWER BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

# OT BOX LARGE

OVERSIZED CASTABLE HOUSING FOR CHAIRSIDE POSITIONING OF CAPS



Ref. 253BLN

## OT BOX LARGE NORMAL & CONNECTORS

### KIT CONTAINS

- 2 OT BOX LARGE BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS



Ref. 253BLM

## OT BOX LARGE MICRO & CONNECTORS

### KIT CONTAINS:

- 2 OT BOX LARGE BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

# OT STRATEGY ATTACHMENTS



Ref. 098SSS

## OT STRATEGY

### CAPS FOR DUPLICATION TECHNIQUE

#### KIT CONTAINS

- 4 CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH)
- 2 CASTABLE STEADY
- 4 RETENTIVE CAPS FOR DUPLICATION TECHNIQUE



Ref. 098CAL

## OT STRATEGY

### CAPS FOR STAINLESS STEEL HOUSING TECHNIQUE

#### KIT CONTAINS

- 4 CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH)
- 2 CASTABLE STEADY
- 4 STAINLESS STEEL HOUSINGS
- 2 PLASTIC POSITIONERS
- 4 RETENTIVE CAPS FOR HOUSING TECHNIQUE

# OT BAR MULTIUSE ATTACHMENTS



Ref. 021OBM

## OT BAR MULTIUSE & CONNECTORS

### KIT CONTAINS

- 2 CASTABLE BARS
- 4 POSITIONER CLIP A
- 4 POSITIONER CLIP B
- 4 CASTABLE BOXES
- 4 RETENTIVE PINK CLIPS
- 4 RETENTIVE YELLOW CLIPS
- 4 RETENTIVE RED CLIPS
- 2 CONNECTORS
- 1 GINGIVAL CONNECTOR

# OT VERTICAL ATTACHMENTS



Ref. 071OBV

## OT VERTICAL

### KIT CONTAINS

- 4 CASTABLE MALES
- 2 CASTABLE STEADY
- 4 RETENTIVE WHITE CLIPS
- 2 RETENTIVE GREEN CLIPS
- 4 CERAMIC BARS
- 4 CASTABLE PARALLELOMETER KEYS + PIN





Ref. 011SPL

## S.P.L. PIVOTS - COPING COVER

### PROMOTIONAL KIT FOR CLINIC

#### KIT CONTAINS

- 24 CASTABLE ASSORTED PIVOTS MICRO SIZE
- 40 CASTABLE ASSORTED PIVOTS NORMAL SIZE
- 12 TITANIUM MICRO PIVOTS
- 12 TITANIUM NORMAL PIVOTS
- 16 NORMAL CAPS (Clear - standard retention - PINK soft retention)
- 16 MICRO CAPS (Clear - standard retention - PINK soft retention)
- 20 PROTECTIVE DISKS
- 2 REGULATING TOOLS (1 normal and 1 micro size)
- 1 BLUE PLASTIC "MULTIUSE" INSERTION TOOL  
(FOR THE INSERTION OF THE CAPS INSIDE THE METAL HOUSINGS)
- 1 OT CAP insertion tool
- 4 NORMAL STAINLESS STEEL HOUSING
- 4 MICRO STAINLESS STEEL HOUSING
- 12 DIRECTIONAL RINGS ( grey 0° - green 7° - red 14°)
- 1 TWEEZERS



Ref. 033PSF

## PIVOT FLEX TiN COATED

TITANIUM PIVOT WITH SWIVEL SPHERE - Normal size  
Adapted for COPING COVER

#### KIT CONTAINS

- 1 TITANIUM PIVOT WITH SWIVEL SPHERE
- 1 STAINLESS STEEL HOUSING
- 2 PINK CAPS Normal Size - Soft retention
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

Available for lengths  
Normal sphere ø 2,5  
Micro sphere ø 1,8



Ref. 036PTN

## TITANIUM PIVOTS KIT NORMAL SIZE

Adapted for COPING COVER

#### KIT CONTAINS

- 5 TITANIUM PIVOTS (Adapted for COPING COVER)

Ref. 036PTM

## TITANIUM PIVOTS KIT MICRO SIZE

Adapted for COPING COVER

#### KIT CONTAINS

- 5 TITANIUM PIVOTS (Adapted for COPING COVER)

STANDARD LENGTH



Ref. 034PRK

## OVERDENTURE POSTS

FOR ROOTS, AVAILABLE IN ONE STANDARD SIZE

#### KIT CONTAINS

- 2 Root posts in Titanium + TiN
- 2 Retentive male attachment Titanium + nylon
- 2 Plastic manual tools
- 2 Stainless steel housings
- 2 Pink protective disks



Sizes available:  
mm. 14 - 12 - 10 - 9 - 7

Ref. 010PSP

## CASTABLE PIVOTS

FOR CANAL IMPRESSIONS

#### KIT CONTAINS

- 20 CASTABLE PIVOTS NORMAL SIZE



Sizes available:  
mm. 10 - 9 - 7

Ref. 012PSM

## CASTABLE PIVOTS

FOR CANAL IMPRESSIONS

#### KIT CONTAINS

- 20 CASTABLE PIVOTS MICRO SIZE



Ref. A01MOG

## MOOSER BUR

FOR PIVOTS OF 7 - 9 - 10mm.



Ref. A03MOB

## MOOSER BUR

FOR PIVOTS OF 12 - 14mm.



Ref. 121FD

OT REVERSE 3 DIAMOND SHAPED BUR FOR ROOTS

**SPHERO FLEX WITH TiN COATED**



**Ref. 109**

**KIT CONTAINS**

- 1 TITANIUM ABUTMENT WITH SWIVEL SPHERE (SPHERE Ø 2.5 mm.)
- 3 CAPS (Yellow - Pink - Clear)
- 1 TITANIUM HOUSING
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

**SPHERO BLOCK NORMO WITH TiN COATED**



**Ref. 002**

**KIT CONTAINS**

- 1 TITANIUM ABUTMENT (SPHERE Ø mm. 2.5)
- 3 CAPS (Yellow - Pink - Clear)
- 1 TITANIUM HOUSING
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

**SPHERO BLOCK MICRO WITH TiN COATED**



**Ref. 003**

**KIT CONTAINS:**

- 1 TITANIUM ABUTMENT (SPHERE Ø mm. 1.8)
- 3 CAPS (Yellow - Pink - Clear)
- 1 TITANIUM HOUSING
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

THE SPHERO FLEX, SPHERO BLOCK ATTACHMENTS ARE AVAILABLE FOR THE MAJOR IMPLANT BRANDS. UPON REQUEST WE CAN CUSTOMIZE THEM FOR ANY KIND OF IMPLANTS. WHEN ORDERING, PLEASE ALWAYS INDICATE: IMPLANT BRAND, IMPLANT DIAMETER AND CUFF HEIGHT. FOR SPHERO FLEX ABUTMENT CUFF AVAILABLE FROM 0,5 TO 7mm.



**AUXILIARY TOOLS SPHERO FLEX - BLOCK**



**Ref. 771CEF**

**KEY FOR SPHERO FLEX AND SPHERO BLOCK (Normal - Micro size)**  
HEX mm. 2.3



**Ref. 760CBM**

**HEX SCREWDRIVER FOR TORQUE CONTROLLER CONTRA ANGLE**



**Ref. 760CRD**

**RATCHET TORQUE CONTROL DEVICE**  
For Sphero block - flex and Ot Equator  
15/35Ncm Strength - Max 50Ncm torque, suggested 25Ncm.

**ANTI-UNSCREWING SYSTEM WITH AN ELASTIC DOWEL APPLICABLE TO**

SPHERO FLEX-BLOCK - OT REVERSE AND DIFFERENT SCREWS. THE ELASTIC DOWEL IS APPLICABLE ONLY UPON REQUEST TO ALL SCREWS WITH A NOT INFERIOR M. 1.8 to screws - Sphero Flex - Sphero Block - OT Reverse attachments, etc.

**Ref. TE**

Extra charge for the dowel application



## AUXILIARY TOOLS SPHERO FLEX - BLOCK



Ref. 100AD

### DIRECTIONAL RINGS FOR SPHERO FLEX AND BLOCK

#### KIT CONTAINS

- 3 COLOURED RINGS WITH INCLINATION OF 0° - 7° - 14°



Ref. 044ABN

### SPHERO FLEX - BLOCK NORMO ANALOG KIT

#### KIT CONTAINS

- 1 SPHERO BLOCK-FLEX NORMO ANALOG
- 3 DIRECTIONAL RINGS (0° - 7° - 14°)
- 1 BLACK NORMO CAP FOR LABORATORY USE



Ref. 044ABM

### SPHERO BLOCK MICRO ANALOG KIT

#### KIT CONTAINS:

- 1 SPHERO BLOCK-MICRO ANALOG
- 3 DIRECTIONAL RINGS (0° - 7° - 14°)
- 1 BLACK MICRO CAP FOR LABORATORY USE

## CUFF HEIGHT MEASURER

THIS TOOL IS USED TO DETERMINE THE CUFF HEIGHT OF IMPLANTS WITH INTERNAL OR EXTERNAL HEXAGON



Ref. 008MBG

### CUFF HEIGHT MEASURER

#### KIT CONTAINS:

- 1 CUFF HEIGHT SLIDER GAUGE
- 1 CUFF HEIGHT FIXED ROD GAUGE
- 1 SILICON RING DISPENSER
- 20 SILICON RINGS



## OT LOCK

LOCKING PIN TITANIUM AND CASTABLE



Ref. 880CLT

### OT LOCK KIT TITANIUM

#### KIT CONTAINS

- 1 LOCKING PIN SET
- 1 HOUSING SHAPER
- 1 CERAMIC PIN



Ref. 880CLR

### OT LOCK KIT TITANIUM ADJUSTABLE

#### KIT CONTAINS

- 1 LOCKING PIN SET
- 1 HOUSING SHAPER
- 1 CERAMIC PIN
- 9 CASTABLE SPACER RINGS

## OT LOCK ACCESSORIES



Ref. 882CG

- 1 CONICAL GUIDE



Ref. 882CAS

- 1 UNLOCKING TOOL

# SPARE PARTS OT CAP SYSTEM NORMAL Ø2.5mm



- Ref. 055SCN** • 4 CASTABLE SINGLE SPHERES NORMAL SIZE
- Ref. 151BAN** • 2 CASTABLE PRESECTIONED BEVELED BARS NORMAL SIZE
- Ref. 150BPN** • 2 CASTABLE PRESECTIONED BARS NORMAL SIZE



- Ref. 192ACN** **CAPS ASSORTMENT KIT OT CAP NORMAL SIZE**  
**KIT CONTAINS**
  - 1 TITANIUM HOUSINGS
  - 1 PROCESSING CAP
  - 1 PROTECTIVE DISK
  - 4 RETENTIVE CAPS (1 GREEN, 1 YELLOW, 1 PINK, 1 CLEAR)



- Ref. 048CAN** • 6 SILVER CAPS NORMAL SIZE (EXTRA-RESILIENT) (ELASTIC AND GUMMY) - 350 g



- Ref. 048CON** • 6 GOLD CAPS NORMAL SIZE (EXTRA-RESILIENT) (SLIGHTLY ELASTIC) - 500g



- Ref. 049PCN** • 6 GREEN CAPS NORMAL SIZE (VERY ELASTIC RETENTION) - 350g



- Ref. 060CRN AY** • 6 YELLOW CAPS NORMAL SIZE (EXTRA SOFT RETENTION) - 500g



- Ref. 040CRN SN** • 6 PINK CAPS NORMAL SIZE (SOFT RETENTION) - 900g



- Ref. 040CRN** • 6 CLEAR CAPS NORMAL SIZE (STANDARD RETENTION) - 1300g



- Ref. 040TCN** • 2 TITAN CAPS NORMAL SIZE



- Ref. 043CLN** • 6 BLACK CAPS NORMAL SIZE (FOR LABORATORY USE ONLY)



- Ref. 049PCNDR8** • 6 ORANGE CAPS NORMAL SIZE (VERY ELASTIC RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 350 gr



- Ref.060CRNAYDR8** • 6 YELLOW CAPS NORMAL SIZE (EXTRA SOFT RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 500 gr



- Ref.040CRNSNDR8** • 6 PINK CAPS NORMAL SIZE (SOFT RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 900 gr



- Ref. 040CRNDR8** • 6 AQUA CAPS NORMAL SIZE (STRONG RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 1300 gr



- Ref. 042AMN** • 2 MONO OT BOX NORMAL SIZE (CASTABLE HOUSINGS FOR CAPS)



- Ref. 041CAN** • 2 STAINLESS STEEL HOUSINGS NORMAL SIZE



- Ref. 041CTN** • 2 TITANIUM HOUSINGS NORMAL SIZE



- Ref. 101PPN** • 2 PLASTIC POSITIONERS NORMAL SIZE



- Ref. 044PPN** • 2 NORMAL SPHERE ANALOGS



- Ref. 100SAN** • 10 ALUMINIUM SPACERS FOR OT BOX BARS NORMAL SIZE

## SPARE PARTS OT CAP SYSTEM MICRO Ø1.8mm



**Ref. 055SCM**

• 4 CASTABLE SPHERES MICRO SIZE

**Ref. 151BAM**

• 2 CASTABLE PRESECTIONED BEVELED BARS MICRO SIZE

**Ref. 150BPM**

• 2 CASTABLE PRESECTIONED BARS MICRO SIZE



**Ref. 192ACM**

### CAPS ASSORTMENT KIT OT CAP MICRO SIZE

#### KIT CONTAINS

- 1 TITANIUM HOUSINGS
- 1 PROCESSING CAP
- 1 PROTECTIVE DISK
- 4 RETENTIVE CAPS (1 GREEN, 1 YELLOW, 1 PINK, 1 CLEAR)



**Ref. 048CAM**

• 6 SILVER CAPS MICRO SIZE (EXTRA-RESILIENT)  
(ELASTIC AND GUMMY) - 200g



**Ref. 048COM**

• 6 GOLD CAPS MICRO SIZE (EXTRA-RESILIENT)  
(SLIGHTLY ELASTIC) - 450g



**Ref. 049PCM**

• 6 GREEN CAPS MICRO SIZE  
(VERY ELASTIC RETENTION) - 200g



**Ref. 060CRM AY**

• 6 YELLOW CAPS MICRO SIZE  
(EXTRA SOFT RETENTION) - 450g



**Ref. 040CRM SN**

• 6 PINK CAPS MICRO SIZE  
(SOFT RETENTION) - 800g



**Ref. 040CRM**

• 6 CLEAR CAPS MICRO SIZE  
(STANDARD RETENTION) - 1100g



**Ref. 040TCM**

• 2 TITAN CAPS MICRO SIZE



**Ref. 043CLM**

• 6 BLACK CAPS MICRO SIZE (FOR LABORATORY USE ONLY)



**Ref. 040CRMDR8**

• 6 AQUA CAPS MICRO SIZE (STRONG RETENTION)  
(UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø1,6 mm) - 1100g



**Ref. 049PCMDR8**

• 6 ORANGE CAPS MICRO SIZE (VERY ELASTIC RETENTION)  
(UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø1,6 mm) - 200g



**Ref. 042AMM**

• 2 MONO OT BOX MICRO SIZE  
(CASTABLE HOUSINGS FOR CAPS)



**Ref. 041CAM**

• 2 STAINLESS STEEL HOUSINGS MICRO SIZE



**Ref. 041CTM**

• 2 TITANIUM HOUSINGS MICRO SIZE



**Ref. 101PPM**

• 2 PLASTIC POSITIONERS MICRO SIZE



**Ref. 044PPM**

• 2 MICRO SPHERE ANALOGS



**Ref. 100SAM**

• 10 ALUMINIUM SPACERS FOR OT BOX BARS MICRO SIZE

## NORMAL-MICRO TOOLS



**Ref. 485IC**

• EXTRACTOR/INSERTER TOOL FOR CAPS - NORMAL AND MICRO SIZES  
(FOR THE INSERTION OF THE CAPS INSIDE THE METAL HOUSINGS)



**Ref. 74AC01**

• KEY FOR PARALLELOMETER NORMAL SIZE



**Ref. 74AC02**

• KEY FOR PARALLELOMETER MICRO SIZE



**Ref. 74AC03**

• TECNO KEY FOR PARALLELOMETER NORMAL AND MICRO SIZES  
(ONLY FOR OT CAP TECNO)



**Ref. 491EC**

• MULTIUSE REMOVAL TOOL FOR RETENTIVE CAPS

## NORMAL-MICRO TOOLS continue



Ref. 080RCN

- REAMER TOOL FOR OT CAPS - NORMAL SIZE



Ref. 080RCM

- REAMER TOOL FOR OT CAPS - MICRO SIZE



Ref. 082ATN

- TOOL FOR TESTING CAP RETENTION - NORMAL SIZE



Ref. 083ATM

- TOOL FOR TESTING CAP RETENTION - MICRO SIZE

## NORMAL - MICRO ACCESSORIES AND Ø2.2mm



Ref. 155CC

- 6 CONNECTORS



Ref. 187CRS

- 2 CONCAVE SPHERES  
AVAILABLE: Ø1.8 - 2.2 - 2.5mm



Ref. 187CRE

- 2 CONCAVE OT EQUATOR



Ref. 189SRS

- 2 SOLID SPHERES  
AVAILABLE: Ø1.8mm



Ref. 100PD

- 10 PROTECTIVE DISKS (internal diameter 1,6mm)



Ref. 100PDR

- 10 PINK PROTECTIVE DISKS (internal diameter 2,9mm)



Ref. 038GCS

- 2 CASTABLE SLIDING BASES



Ref. 044PP22

- 2 MICRO ANALOGS Ø2.2mm

## OT CAP ASSORTMENT KIT

### KIT NORMAL

Ref. 074KAN

### KIT MICRO

Ref. 074KAM

### KIT CONTAINS

- 8 CLEAR CAPS - STANDARD RETENTION
- 8 PINK CAPS - SOFT RETENTION
- 6 YELLOW CAPS - EXTRA SOFT RETENTION
- 6 GREEN CAPS - VERY ELASTIC RETENTION
- 6 GOLD CAPS - (EXTRA RESILIENT) SLIGHTLY ELASTIC
- 6 SILVER CAPS - (EXTRA RESILIENT) ELASTIC AND GUMMY
- 4 LIGHT GREEN CAPS - STANDARD (FOR WORN OUT SPHERE)
- 4 ORANGE CAPS - ELASTIC (FOR WORN OUT SPHERE)
- 4 BLACK CAPS - LABORATORY USE ONLY
- 2 TITAN CAPS

### ASSORTMENT KIT NORMAL & MICRO

#### KIT CONTAINS

- 4+4 CLEAR CAPS - STANDARD RETENTION
- 4+4 PINK CAPS - SOFT RETENTION
- 2+2 YELLOW CAPS - EXTRA SOFT RETENTION
- 3+3 GREEN CAPS - VERY ELASTIC RETENTION
- 3+3 GOLD CAPS - (EXTRA RESILIENT) SLIGHTLY ELASTIC
- 3+3 SILVER CAPS - (EXTRA RESILIENT) ELASTIC AND GUMMY
- 2+2 LIGHT GREEN CAPS - STANDARD (FOR WORN OUT SPHERE)
- 2+2 ORANGE CAPS - ELASTIC (FOR WORN OUT SPHERE)
- 2+2 BLACK CAPS - LABORATORY USE ONLY
- 1+1 TITAN CAPS

### OT CAP ASSORTED RETENTIVE CAPS

#### NORMAL KIT

- 6 GREEN CAPS - VERY ELASTIC RETENTION
- 6 YELLOW CAPS - EXTRA SOFT RETENTION
- 6 PINK CAPS - SOFT RETENTION
- 6 CLEAR CAPS - STANDARD RETENTION

### OT CAP ASSORTED RETENTIVE CAPS

#### MICRO KIT

- 6 GREEN CAPS - VERY ELASTIC RETENTION
- 6 YELLOW CAPS - EXTRA SOFT RETENTION
- 6 PINK CAPS - SOFT RETENTION
- 6 CLEAR CAPS - STANDARD RETENTION



# OT STRATEGY SPARE PARTS



- Ref. 154PCS** • 4 OT STRATEGY STANDARD MALES
- Ref. 054PLS** • 4 OT STRATEGY HIGH MALES
- Ref. 154STY** • 2 CASTABLE STEADY
- Ref. 047CSD AY** • 4 STRATEGY YELLOW CAPS FOR DUPLICATION TECHNIQUE  
(EXTRA SOFT RETENTION) - 500gr
- Ref. 047CSD SN** • 4 STRATEGY PINK CAPS FOR DUPLICATION TECHNIQUE  
(SOFT RETENTION) - 900gr
- Ref. 047CSD** • 4 STRATEGY CLEAR CAPS FOR DUPLICATION TECHNIQUE  
(STANDARD RETENTION) - 1300gr
- Ref. 047CLS** • 4 STRATEGY BLACK CAPS FOR DUPLICATION TECHNIQUE  
(FOR LABORATORY USE ONLY)



- Ref. 047CSDDR8** • 4 AQUA STRATEGY CAPS FOR DUPLICATION TECHNIQUE  
(UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE)



- Ref. 045CSM AY** • 4 STRATEGY YELLOW CAPS FOR S/STEEL HOUSING  
(EXTRA SOFT RETENTION) - 500gr  
(TO BE USED ONLY WITH HOUSINGS)
- Ref. 045CSM SN** • 4 STRATEGY PINK CAPS FOR S/STEEL HOUSING  
(SOFT RETENTION) - 900gr  
(TO BE USED ONLY WITH HOUSINGS)
- Ref. 045CSM** • 4 STRATEGY CLEAR CAPS FOR S/STEEL HOUSING  
(STANDARD RETENTION) - 1300gr  
(TO BE USED ONLY WITH HOUSINGS)
- Ref. 045CLS** • 4 STRATEGY BLACK CAPS FOR S/STEEL HOUSING  
(FOR LABORATORY USE ONLY)



- Ref. 045CSMDR8** • 4 AQUA STRATEGY CAPS FOR S/STEEL HOUSING  
(UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE)



- Ref. 046CSM** • 2 STRATEGY STAINLESS STEEL HOUSINGS



- Ref. 101OTS** • 2 STRATEGY PLASTIC POSITIONERS  
(FOR SOLDERING STAINLESS STEEL HOUSINGS)

## OT STRATEGY ASSORTMENT CAP KIT FOR DUPLICATION TECHNIQUE



- Ref. 047ACS** **KIT CONTAINS**
  - 4 YELLOW CAPS - EXTRA SOFT RETENTION
  - 4 PINK CAPS - SOFT RETENTION
  - 4 CLEAR CAPS - STANDARD RETENTION

## OT STRATEGY ASSORTMENT CAP KIT FOR S/STEEL HOUSINGS



- Ref. 045ACS** **KIT CONTAINS**
  - 4 YELLOW CAPS - EXTRA SOFT RETENTION
  - 4 PINK CAPS - SOFT RETENTION
  - 4 CLEAR CAPS - STANDARD RETENTION



- Ref. 044PPS** • 2 STRATEGY ANALOGS



- Ref. 486ICS** • INSERTER/EXTRACTOR TOOL FOR STRATEGY CAPS



- Ref. 75AC04** • PARALLELOMETER KEY FOR STRATEGY CAPS




- Ref. 081RCS** • REAMER TOOL FOR STRATEGY CAPS

## OT BAR MULTIUSE SPARE PARTS



|                     |   |
|---------------------|---|
| <b>Ref. 022OBB</b>  | • 2 CASTABLE BARS OT BAR MULTIUSE                   |
| <b>Ref. 025CPB</b>  | • 4 CASTABLE BOXES OT BAR MULTIUSE                  |
| <b>Ref. 023CPA</b>  | • 4 POSITIONER CLIPS A                              |
| <b>Ref. 024CPB</b>  | • 4 POSITIONER CLIPS B                              |
| <b>Ref. 027CRG</b>  | • 4 YELLOW CLIPS (MEDIUM RETENTION 1,8Kg)           |
| <b>Ref. 026CRR</b>  | • 4 PINK CLIPS (SOFT RETENTION 1,2Kg)               |
| <b>Ref. 027CRS</b>  | • 4 RED CLIPS (EXTRA SOFT RETENTION 0,6Kg)          |
| <b>Ref. 022RGO</b>  | • 2 CASTABLE CONNECTORS                             |
| <b>Ref. 429IOBM</b> | • INSERTER/EXTRACTOR TOOL FOR OT BAR MULTIUSE CLIPS |
| <b>Ref. 028OCP</b>  | • PARALLELOMETER KEY FOR OT BAR MULTIUSE            |

## OT VERTICAL SPARE PARTS



|                    |  |
|--------------------|--|
| <b>Ref. 072MCV</b> | • 2 CASTABLE MALES<br>• 2 PIN KEYS       |
| <b>Ref. 072SCV</b> | • 2 CASTABLE STEADY                      |
| <b>Ref. 072BCV</b> | • 2 CERAMIC BARS                         |
| <b>Ref. 072CBV</b> | • 4 WHITE CLIPS (STANDARD RETENTION)     |
| <b>Ref. 072CVV</b> | • 4 GREEN CLIPS (SOFT RETENTION)         |
| <b>Ref. 472ICV</b> | • 1 INSERTION TOOL FOR OT VERTICAL CLIPS |

## OT REVERSE 3 SPARE PARTS



|                    |   |
|--------------------|---|
| <b>Ref. 120MR2</b> | • 2 RETENTIVE MALES in TITANIUM + NYLON (Color Clear/ STANDARD RETENTION) |
| <b>Ref. 041CAM</b> | • 2 S/STEEL HOUSINGS MICRO SIZE   |
| <b>Ref. 121AOR</b> | • 2 BRASS ANALOGS FOR ROOTS   |
| <b>Ref. 121AAR</b> | • 2 INOX IMPLANT ANALOGS  |
| <b>Ref. 121MT</b>  | • 2 PLASTIC TRANSFERS FOR IMPRESSION                                      |





**Ref. 121AM**

- 2 PLASTIC MANUAL TOOLS

**Ref. 100PDR**

- 10 PINK PROTECTIVE DISKS

**Ref. 121PR**

- 2 PLASTIC POSITIONERS FOR DUPLICATION

**Ref. 124IMR**

- METAL INSERTING TOOL FOR OT REVERSE/3 MALE

**Ref. 72CBR**

- S/STEEL PARALLELOMETER KEY FOR CASTABLE BAR

**Ref. 72CFR**

- BRASS PARALLELOMETER KEY FOR TITANIUM FEMALE

## OT LOCK SPARE PARTS



**Ref. 882CS**

- 1 SLIDING KEY  
*(ONLY FOR TITANIUM - CASTABLE LOCKING PIN)*



**Ref. 882AF**

- 1 FLEXIBLE RING



**Ref. 882CB**

- 1 LOCKING KEY



**Ref. 882CPT**

- 1 TITANIUM PREFABRICATED HOUSING  
*(ONLY FOR TITANIUM - TITANIUM ADJUSTABLE LOCKING PIN)*



**Ref. 882PPC**

- 1 CERAMIC PIN



**Ref. 882CLSCT**

- 1 ADJUSTABLE SLIDING KEY  
*(ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)*



**Ref. 882DCC**

- 9 CASTABLE SPACER RINGS  
*(ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)*

## SPHERE REPAIR KITS



AVAILABLE Ø  
1.8 - 2.2 - 2.5

**Ref. 087CRS**

### CONCAVE REPAIR SPHERE

#### KIT CONTAINS

- 2 CONCAVE SPHERES IN TITANIUM - TiN COATED
- 2 PINK CAPS SOFT RETENTION
- 1 INSERTING TOOL
- 1 GAUGE AND STRIP HOLDER



**Ref. 087CRE**

### CONCAVE REPAIR OT EQUATOR

#### KIT CONTAINS

- 2 CONCAVE OT EQUATOR IN TITANIUM - TiN COATED
- 2 PINK CAPS SOFT RETENTION
- 1 INSERTING TOOL
- 1 GAUGE AND STRIP HOLDER



**Ref. 089SRS**

### SOLID REPAIR SPHERE

#### KIT CONTAINS

- 2 SOLID SPHERES 1.8mm IN TITANIUM - TiN COATED
- 2 PINK CAPS SOFT RETENTION
- 2 PROTECTIVE DISKS
- 1 KEY

## TITANIUM ATTACHMENTS FOR OVERDENTURE OT REVERSE SYSTEM - OT KLASSE:

**OT REVERSE SPARE PARTS AND OT KLASSE ARE IN STOCK  
FOR DETAILED INFORMATION PLEASE E-MAIL US OR CALL OUR OFFICES**

## ADJUSTABLE DENTURE BOXES

IT IS POSSIBLE TO SET UP SIX DIFFERENT INTERIOR SIZES



Ref. 0P6

### KIT CONTAINS


- 6 ASSORTED BOXES  
COLOURS:  
BLUE - LIGHT BLUE - WHITE



Ref. 100AS

- 30 SILICONE RINGS

## SMALL PARALLELOMETER



Ref. 00PB

- SMALL PARALLELOMETER DEVICE  
WITH MODEL HOLDER BASE  
(FOR COURSES, ETC.)

## OT CEM



Ref. 0C

### KIT CONTAINS

- 1 TWO-COMPONENT SYRINGE 5ml
- 10 AUTOMIX TIPS

## BROKEN SCREW EXTRACTOR KIT

BROKEN SCREW EXTRACTOR KITS ARE MANUFACTURED FOR MOST IMPLANT BRANDS



Ref. 680

### KIT CONTAINS

- 1 CLAW REAMER BUR
- 1 REVERSE CUTTING BUR
- 1 MANUAL EXTRACTOR
- 1 CENTERING DEVICE

Ref. 680FS

- 1 REVERSE CUTTING BUR

Ref. 680FA

- 1 CLAW REAMER BUR

# RHEIN83

# RHEIN83

## DIRECT SERVICES

A telephone service is available during office hours,  
it furnishes direct answers for technical questions.  
You can address to Rhein83 for knowing the address of foreign distributors

**Tel. 051 244510 - 051 244396 - Fax 051 245238**

Rhein83 website  
**[www.rhein83.com](http://www.rhein83.com)**

General informations  
**[info@rhein83.com](mailto:info@rhein83.com)**

Sales department  
**[angela.internationaltrade@rhein83.it](mailto:angela.internationaltrade@rhein83.it)**

Technical support  
**[gianni.storni@rhein83.it](mailto:gianni.storni@rhein83.it)**

Marketing department  
**[marketing@rhein83.it](mailto:marketing@rhein83.it)**





Via ZAGO, 10/ABC  
40128 - BOLOGNA

Tel. (+39) 051 244510 - (+39) 051 244396  
Fax (+39) 051 245238

<http://www.rhein83.com>  
e-mail: [info@rhein83.com](mailto:info@rhein83.com)

Distributor



UNI EN ISO 9001:2008  
UNI CEI EN ISO 13485:2012